

CT Dedicated Renal Protocol with Pelvis

GE- 6011

Maximum CTDI 75

Indication: Renal masses, nephrectomy follow-up

PT Prep: 1,000 ml of water **30 minutes** prior to getting on the table
 IV contrast – Yes (follow IV contrast administration guidelines)

Series 1: Scouts AP & LAT – Supine “O” at Xiphoid Process S20 to I550

Series 2: Unenhanced – Scan from Top of Kidneys to Iliac Crest

	750 HD (128) CT2	Optima 660 (32) OVIC	Optima 660 CT1	VCT 64 CT3
Noise Level	16.40	12.00	11.60	11.60
Interval	2.5 mm	2.5 mm	2.5 mm	2.5 mm
Axial/Helical Thickness	2.5 mm	2.5 mm	2.5 mm	2.5 mm
Pitch	0.984:1	0.984:1	0.984:1	0.984:1
Speed mm/rotation	39.37	39.37	39.37	39.37
Detector Rows				
Detector Configuration				
Beam Collimation	40 mm	40 mm	40 mm	40mm
KV/mA	Auto mA – if large pt. use manual & maximize mA	Auto mA – if large pt. use manual & maximize mA	Auto mA – if large pt. use manual & maximize mA	Auto mA – if large pt. use manual & maximize mA
Scan Type	Helical Full 0.8 sec	Helical Full 0.8 sec	Helical Full 0.8 sec	Helical Full 0.8 sec

Series 3: Corticomedullary Phase: IV contrast per MCR contrast dosing guidelines @ 3cc/sec.

Delay: 30 seconds – Scan from top of kidneys through iliac crest

	750 HD (128) CT2	Optima 660 (32) OVIC	Optima 660 CT1	VCT 64 CT3
Noise Level	15.00	12.00	11.60	11.60
Interval	2.5 mm	2.5 mm	2.5 mm	2.5 mm
Axial/Helical Thickness	2.5 mm	2.5 mm	2.5 mm	2.5 mm
Pitch	0.984:1	1.375:1	0.984:1	0.984:1
Speed mm/rotation	39.37	39.37	39.37	39.37
Detector Rows				
Detector Configuration				
Beam Collimation	40 mm	40 mm	40 mm	40 mm

Original Date: 4-8-04

Approved by: Dr. Songmen, MCR _____

Revised Date: 11/9/04, 1/11/10, 12/8/10, 04/17/13, 12/15/2015, 2/28/18, 07/27/21, 1/30/24

Page 1

CT Dedicated Renal Protocol with Pelvis

GE- 6011

Maximum CTDI 75

KV/mA	Auto mA – if large pt. use manual & maximize mA	Auto mA – if large pt. use manual & maximize mA	Auto mA – if large pt. use manual & maximize mA	Auto mA – if large pt. use manual & maximize mA
Scan Type	Helical Full 0.8 sec	Helical Full 0.8 sec	Helical Full 0.8 sec	Helical Full 0.8 sec

Series 4: Nephrographic Phase:

Delay: 80 seconds – Scan from top of kidneys through symphysis pubis

	750 HD (128) CT2	Optima 660 (32) OVIC	Optima 660 CT1	VCT 64 CT3
Noise Level	15.00	12.00	11.60	11.60
Interval	2.5 mm	2.5 mm	2.5 mm	2.5 mm
Axial/Helical Thickness	2.5 mm	2.5 mm	2.5 mm	2.5 mm
Pitch	0.984:1	1.375:1	0.984:1	0.984:1
Speed mm/rotation	39.37	39.37	39.37	39.37
Detector Rows				
Detector Configuration				
Beam Collimation	40 mm	40 mm	40 mm	40 mm
KV/mA	Auto mA – if large pt. use manual & maximize mA	Auto mA – if large pt. use manual & maximize mA	Auto mA – if large pt. use manual & maximize mA	Auto mA – if large pt. use manual & maximize mA
Scan Type	Helical Full 0.8 sec	Helical Full 0.8 sec	Helical Full 0.8 sec	Helical Full 0.8 sec

Series 5: Excretory Phase

Delay: 240 seconds – Scan from top of kidney through iliac crest

	750 HD (128) CT2	Optima 660 (32) OVIC	Optima 660 CT1	VCT 64 CT3
Noise Level	15.00	12.00	11.60	11.60
Interval	2.5 mm	2.5 mm	2.5 mm	2.5 mm
Axial/Helical Thickness	2.5 mm	2.5 mm	2.5 mm	2.5 mm
Pitch	0.984:1	1.375:1	0.984:1	0.984:1
Speed mm/rotation	39.37	27.50	39.37	39.37
Detector Rows				
Detector Configuration				
Beam Collimation	40 mm	40 mm	40 mm	40 mm
KV/mA	Auto mA – if large pt. use manual & maximize mA	Auto mA – if large pt. use manual & maximize mA	Auto mA – if large pt. use manual & maximize mA	Auto mA – if large pt. use manual & maximize mA
Scan Type	Helical Full 0.8 sec	Helical Full 0.8sec	Helical Full 0.8 sec	Helical Full 0.8 sec

Original Date: 4-8-04

Approved by: Dr. Songmen, MCR _____

Revised Date: 11/9/04, 1/11/10, 12/8/10, 04/17/13, 12/15/2015, 2/28/18, 07/27/21, 1/30/24

Page 2

CT Dedicated Renal Protocol with Pelvis

GE- 6011

Maximum CTDI 75

--	--	--	--	--

- Networking/ PACs:**
- Send scouts
 - Send series 2,3,4,5 Standard Soft Tissue Algorithm
 - Recon and send Venous series only in Bone Algorithm
 - Recon and send series 3 in Lung Algorithm
 - Send MPR's all series
 - Record DLP in PACS comments