

Job Aid: Ultrasound Transvaginal Protocol		
Division: SLH,SNGH	Manual:	Locations: Radiology Ultrasound
Original Date: 8/1/12	Section:	
Owner: Imaging Operations		
Approved By: Richard Thomas, MD, Department Chief		Revision Date: 11/20/13, 3/04/14, 1/20/17, 2/06/2020

Organs to be demonstrated: Cervix, Uterus, Endometrium, Adnexa, and Ovaries.

1. Sagittal Uterus

- 2 images cervix
- 1 image midline uterus
-w/wo calipers if not used at pelvic protocol
- 1 image endometrium (w/wo calipers)
- 2 images RT and LT uterus

2. Transverse Uterus

- 2 images cervix
- 3 images uterus including body and fundus
-w/wo calipers if not used at pelvic protocol

3. SAG/TRV RT Ovary

- 2 images SAG including measurement
- 2 images TRV including measurement
color and pulsed doppler assessment

4. SAG/TRV LT Ovary

- 2 images SAG including meas.
- 2 images TRV including meas.
color and pulsed doppler assessment

5. CUL-DE-SAC

- 2 images SAG cul-de-sac

Imaging Notes:

Do not need to measure all fibroids...Measure largest 1-3 fibroids or any designated For treatment. Show any submucosal and pedunculated fibroids.

When measuring fibroids and/or masses use Dual screen or make sure Sag and Trv Images are consecutive. Measure uterus from cervix to fundus unless there is a Request for uterine volume- measure internal os of cervix to fundus.

***IUD-Please do cine loop in transverse showing correct placement.**

Related Documents: